

Collagen Mass

Product Description

The Collagen Mass is a dispersion of bovine Collagen which is diluted in HCl. It has a viscose, jellylike consistency.

The Collagen Mass is delivered cooled and non-sterile.

Applications

The Collagen Mass can be used for individual research applications especially for coating of implantable devices.



SPECIFICATIONS

Water	85 – 90 %
Collagen	10 – 15 %
pH-value	2,5 – 3

Storage

The Collagen Mass should be stored between +2°C and +8°C in closed packaging. For longer storage it should be stored at -20°C.

Shelf life: 24 months

Collagen Mass is intended for research use only. It is neither intended for human nor animal diagnostic, therapeutic use or any other clinical uses!

All data and recommendations correspond to the present state of our knowledge; they are published without engagement. We reserve the right to make alterations and additions in line with technical developments without prior notice. The customer is obliged to check whether our products meet with his own technical requirements. We shall be glad to answer any queries.

Viscofan BioEngineering

A Business Unit of Naturin Viscofan GmbH
 Badeniastraße 13
 69469 Weinheim
 Germany

Tel.: +49 (0)6201 86-358
 Fax: +49 (0) 6201 86-226
 Email: sales@bio.viscofan.com
www.viscofan-bioengineering.com