

## Collagen Mass

### Product Description

The Collagen Mass is a suspension of fibrous bovine collagen type I which is diluted in HCl. It has a viscose, jellylike consistency.

The Collagen Mass is delivered cooled and non-sterile.

### Applications

The Collagen Mass can be used for individual research applications especially for material coating.



### SPECIFICATIONS

Water	85 – 90 %
Collagen	10 – 15 %
pH-value	2,2 – 3,0

### Storage

The Collagen Mass should be stored between +2°C and +8°C in closed packaging. For longer storage it should be stored at -20°C.

Shelf life: 12 months

*Collagen Mass is intended for research use only. It is neither intended for human nor animal diagnostic, therapeutic use or any other clinical uses!*

All data and recommendations correspond to the present state of our knowledge; they are published without engagement. We reserve the right to make alterations and additions in line with technical developments without prior notice. The customer is obliged to check whether our products meet with his own technical requirements. We shall be glad to answer any queries.

### Viscofan BioEngineering

A Business Unit of Naturin Viscofan GmbH  
 Badeniastraße 13  
 69469 Weinheim  
 Germany

Tel.: +49 (0)6201 86-358  
 Fax: +49 (0) 6201 86-226  
 Email: [sales@bio.viscofan.com](mailto:sales@bio.viscofan.com)  
[www.viscofan-bioengineering.com](http://www.viscofan-bioengineering.com)